

# Labotec Glassware Drying Oven



## Features:

This simple, robust unit is ideally suited for the rapid and efficient drying of laboratory glassware.

- Grade 304 stainless steel body construction
- Reliable incolloy heating elements
- Unit provided with 2 chrome plated wire shelves
- Toughened glass sliding doors on both sides for ease of access
- Variable heat control via energy regulator

## Technical Data – Circulators

Model Code	Electrical Connection	Capacity (L)	External Dimensions (W x D x H mm)	Internal Dimensions (W x D x H mm)	Approx. Weight (Kg)
310	230V 50hz	180	1035 x 510 x 518	965 x 410 x 510	27

