

LABOTEC OVERHEAD STIRRER

EcoStir
LABOTEC

The Labotec EcoStir overhead stirrers are suitable for most types of laboratory applications where controlled stirring is required.

Features

- Maximum volume 20 litres
- Easy to set digital control with rotating knob
- Plate size 200mm x 300mm
- Solid state electronic speed control
- RPM 100 - 2500
- Adjustable height and angle via clamping bar

Technical Data

Model No.	720
Voltage (V/hz)	100-240/50
Power (W)	60
External dimension (mm)	
Width	200
Depth	380
Height	860
*Glass flask not included	
Shaft Diameter	8mm
Chuck Diameter	up to 10mm
Speed range	100-2500rpm
Torque	up to 40Ncm
Viscosity	0-10000mPas



Model 720